



Donated Chemical Probe

*TRPA1 Inhibitor*  
*Probe BAY-390*




**A. Laux-Biehlmann,**  
**S.Mesch, H. Miyataka**  
**Ondozaal** on behalf of the team

**September, 2020**





## Background & Rationale

**TRPA1, a damage sensor activated by a wide range of endogenous mediators**




**TRPA1 activation induces pain, central sensitization, neurogenic inflammation and mast cell-driven inflammation**

TRPA1 in human endometriotic lesions



# Commercially available TRPA1 antagonists

A-079 donated as chemical probe


Origin	Code	Structure	hTRPA1	r   m TRPA1	References
Abbvie	A-079		67 nM (FLIPR Ca <sup>2+</sup> assay)	289   298 nM (FLIPR Ca <sup>2+</sup> assay)	PMID: 22319196 PMID: 21402443

## Reference of TRPA1 probe:

<https://www.sgc-ffm.uni-frankfurt.de/chemProbes#!specificprobeoverview/A-079>

<https://www.sgc-ffm.uni-frankfurt.de/APP/connector/0/146/url/A-079+Antagonist+for+TRPA1.pdf>

## Technical Profile



Absolute Stereochemistry: *R,R* by VCD measurement

### Pharmacology

hTRPA1 FLIPR IC <sub>50</sub> [nM]	16
hTRPA1 Ephys IC <sub>50</sub> [nM]	82
rTRPA1 FLIPR IC <sub>50</sub> [nM]	63
rDRG EPhys IC <sub>50</sub> [nM]	35
m   gp   dog   monkey TRPA1 IC <sub>50</sub> [nM]	73   68   81   19

### Selectivity/ Safety

TRP selectivity IC <sub>50</sub> [μM]	> 25   > 25   >25
hTRPV1   V4   TRPC3   C5   C6	5.6   > 25
hCa <sub>v</sub> 3.2   TASK-3 IC <sub>50</sub> [μM]	> 25   > 30
hERG IC <sub>50</sub> [μM]	> 10
Panlabs Lead Profiling [μM]   Cardiac Profiler [%] inhibition   transactivation Erα @10 uM	Ki hPR-B: 4.0; Ki hDAT: 0.9   no findings   EC <sub>50</sub> ERα: 2.1
Ames MPF	negative

Cardiac profiler: hNav1.5; hCav1.2; hKir2.1; hKvLQT1; hKv4.3

### Molecular Properties

MW <sub>corr</sub> [g/mol]	277
TPSA [Å <sup>2</sup> ]	32
logD 7.5 (calc/exp)	3.3 / 3.3

### PhysChem

Sol <sub>Flask</sub> pH 6.5, 4, 2 [mg/L]	n.d. (oil)
Plasmastab (4h) hum/rat [%]	n.d.
Stability (24h) pH1/7/10 [%]	95/100/100


### Pharmacokinetics

	CL <sub>bl</sub> [L/h/kg]		Fmax [%]			
in vitro PK (hHep/rHep)	0.38	2.4	71	42		
	V <sub>ss</sub> [L/kg]	CL <sub>bl</sub> [L/h/kg]	t <sub>1/2</sub> [h] iv	F [%]		
in vivo PK	2.1	1.4	2.0	38		
	Hum [%]	Rat [%]	Mouse [%]	Williams E [%]		
Protein binding fu	11.5	6.0	8.1	57		
	A-B [nm/s]		B-A [nm/s]		Efflux	
CaCo	293				0.63	
	1A2	2C8	2C9	2D6	3A4	TDI
CYP Inhibition [μM]	>20	>20	>20	>20	>20	no
	PXR		NOEL 1A2 [μg/L]		NOEL 3A4 [μg/L]	
CYP Induction	YELLOW		10000		30000	

m: mouse; gp: guinea pig; DRG: dorsal root ganglion



## Technical Profile



### Pharmacology

hTRPA1 FLIPR (WUP) IC <sub>50</sub> [nM]	> 25000
hTRPA1 Ephys (Evo) IC <sub>50</sub> [nM]	n.d.
rTRPA1 FLIPR (WUP) IC <sub>50</sub> [nM]	> 25000
rDRG EPhys (Evo) IC <sub>50</sub> [nM]	n.,.d.

### Selectivity/ Safety

TRP selectivity IC <sub>50</sub> [μM] hTRPV1   V4   TRPC5   C6	>25   >25   10   >30
hCa <sub>v</sub> 3.2   TASK-3 IC <sub>50</sub> [μM]	n.d.
hERG IC <sub>50</sub> [μM]	> 10
Panlabs Lead Profiling [μM]   Cardiac Profiler [%] inhibition	n.d.   all > 10
Ames MPF	n.d.

Cardiac profiler: hNav1.5; hCav1.2; hKir2.1; hKvLQT1; hKv4.3

### Molecular Properties

MW <sub>corr</sub> [g/mol]	307
TPSA [Å <sup>2</sup> ]	41
logD 7.5 (calc/exp)	3.2 / 2.9

### PhysChem


Sol <sub>Flask</sub> pH 6.5, 4, 2 [mg/L]	n.d. (oil)
Plasmastab (4h) hum/rat [%]	n.d.
Stability (24h) pH1/7/10 [%]	n.d.

### Pharmacokinetics

11.	CL <sub>bl</sub> [L/h/kg]		Fmax [%]			
in vitro PK (hHep/rHep)	0.85	3.5	36	17		
	V <sub>ss</sub> [L/kg]	CL <sub>bl</sub> [L/h/kg]	t <sub>1/2</sub> [h] iv	F [%]		
in vivo PK	n.d.	n.d.	n.d.	n.d.		
	Hum [%]	Rat [%]	Mouse [%]	Williams E [%]		
Protein binding fu	n.d.	n.d.	n.d.	n.d.		
	A-B [nm/s]		B-A [nm/s]		Efflux	
CaCo	375				0.55	
	1A2	2C8	2C9	2D6	3A4	TDI
CYP Inhibition [μM]	> 20	> 20	> 20	> 20	> 20	
	PXR		NOEL 1A2 [μg/L]		NOEL 3A4 [μg/L]	
CYP Induction	RED		30000		2610	


# BAY-390 shows efficacy in inflammatory pain and neurogenic inflammation models

## Reduction of mechanical hyperalgesia and inflammation in mouse inflammatory pain models



**Inflammatory Pain and Peripheral Inflammation**


## Reduction of neurogenic inflammation (plasma extravasation) in a rat neurogenic inflammation model



**Neurogenic Inflammation**


# BAY-390 shows efficacy in neuropathic pain and interstitial cystitis models

## Reduction of mechanical hyperalgesia in rodent neuropathic pain model



**Neuropathic Pain**

## Reduction of visceral pain in a rat cyclophosphamide induced cystitis model



**IC/BPS**

Is a potent, across species active and brain penetrating TRPA1 antagonist

Probe criteria	
Inhibitor/agonist potency: goal is < 100 nM (IC <sub>50</sub> , Kd)	Potency as inhibitor of TRPA1 ion channel demonstrated in FLIPR Ca <sup>2+</sup> assay (CHO or HEK): human TRPA1 IC <sub>50</sub> = 16 nM, dog IC <sub>50</sub> = 19 nM, rat IC <sub>50</sub> = 63 nM; Equipotent on rat, mouse, guinea pig and monkey TRPA1
Selectivity within target family: goal is >30-fold	Surpasses criteria; Selectivity against family members was tested, all >30fold
Selectivity outside target family: describe the off-targets (which may include both binding and functional data)	Surpasses criteria; Selectivity in Eurofins Lead Profiling Screen, GPCR Profiling Screen and Bayer Kinase Panel was performed (only relevant Ki: hDAT 0.9 μM; hPR-B: 4 μM    Erα EC <sub>50</sub> : 2.1 μM )
On target cell activity for cell-based targets: goal is < 1 micromolar IC <sub>50</sub> /EC <sub>50</sub>	Surpasses criteria; Ephys human TRPA1 IC <sub>50</sub> at 82 nM (Patchliner, CHO) Efficacy in several in vitro and in vivo mechanistic models demonstrated; efficacy in PD pain inflammation model
On target cell activity for secreted targets: appropriate alternative such as mouse model or other mechanistic biological assay, e.g., explant culture	Surpasses criteria; Suitable pharmacokinetic profile for in vivo studies in rodents
Neg ctrl: in vitro potency – > 100 times less; Cell activity – >100 times less potent than the probe	Surpasses criteria; Structure related compound BAY-9897 with high micromolar TRPA1 activity (human TRPA1 IC <sub>50</sub> >25 μM (CHO, Ca <sup>2+</sup> ))

**We ask for acceptance of TRPA1 inhibitor BAY-390 as chemical probe,  
accompanied by BAY-9897 as negative control**





# Acknowledgement



## **Bayer:**

Daniel Basting  
Stefan Bäuerle  
Marina Bairlein  
Celine Bordin  
Ulrich Bothe  
Hans Briem  
Bernd Elger  
Christa Hegele-Hartung  
Herbert Himmel  
Andreas Sutter  
Alexis Laux-Biehlmann  
Hideki Miyatake Ondoababal  
Michaele Peters  
Stefanie Mesch  
Andrea Rotgeri \*  
Antje Rottmann  
Andreas Steinmeyer  
Thomas Zollner

## **Evotec:**

Susan Boyce  
Anne-Marie Coelho  
Stuart Flanagan  
Stephen Hess  
Schanila Nawaz  
Daryl Walter

\* Now at Nuvisan




*Thank you*



Questions?





## BAY-390

### Experimental Results

Cat #	Assay Name	Batch*	Spec.	Rep.	Conc.	% Inh.
<b>Compound: CHH039-2017, PT #: 1207556</b>						
107000	Aldose Reductase	401131	rat	2	10 µM	3
107710	ATPase, Na <sup>+</sup> /K <sup>+</sup> , Heart, Pig	401160	pig	2	10 µM	3
112020	Carbonic Anhydrase II	401132	hum	2	10 µM	7
104010	Cholinesterase, Acetyl, ACES	401128	hum	2	10 µM	-2
116020	Cyclooxygenase COX-1	401157	hum	2	10 µM	8
118010	Cyclooxygenase COX-2	401158	hum	2	10 µM	23
124010	HMG-CoA Reductase	401134	hum	2	10 µM	8
132000	Leukotriene LTC <sub>4</sub> Synthase	401133	gp	2	10 µM	-15
199017	Lipoxygenase 15-LO	401198	hum	2	10 µM	23
140010	Monoamine Oxidase MAO-A	401130	hum	2	10 µM	3
140120	Monoamine Oxidase MAO-B	401159	hum	2	10 µM	7
142000	Nitric Oxide Synthase, Neuronal (nNOS)	401135	rat	2	10 µM	-3
199010	Nitric Oxide Synthetase, Inducible (iNOS)	401197	mouse	2	10 µM	5
107300	Peptidase, Angiotensin Converting Enzyme	401129	rabbit	2	10 µM	3
152000	Phosphodiesterase PDE3	401193	hum	2	10 µM	3
154000	Phosphodiesterase PDE4	401194	hum	2	10 µM	-11
156000	Phosphodiesterase PDE5	401393	hum	2	10 µM	7
194020	Thromboxane Synthase	401136	hum	2	10 µM	6
200510	Adenosine A <sub>1</sub>	401253	hum	2	10 µM	18
200610	Adenosine A <sub>2A</sub>	401253	hum	2	10 µM	3
200720	Adenosine A <sub>3</sub>	401176	hum	2	10 µM	4
203100	Adrenergic α <sub>1A</sub>	401240	rat	2	10 µM	18
203630	Adrenergic α <sub>2A</sub>	401239	hum	2	10 µM	7
203710	Adrenergic α <sub>2B</sub>	401142	hum	2	10 µM	2
203810	Adrenergic α <sub>2C</sub>	401211	hum	2	10 µM	-1
204010	Adrenergic β <sub>1</sub>	401230	hum	2	10 µM	3
204110	Adrenergic β <sub>2</sub>	401239	hum	2	10 µM	11
204200	Adrenergic β <sub>3</sub>	401359	hum	2	10 µM	1
206000	Androgen (Testosterone)	401322	hum	2	10 µM	47
210030	Angiotensin AT <sub>1</sub>	401212	hum	2	10 µM	4
210120	Angiotensin AT <sub>2</sub>	401213	hum	2	10 µM	1
212510	Bradykinin B <sub>1</sub>	401177	hum	2	10 µM	9
212620	Bradykinin B <sub>2</sub>	401170	hum	2	10 µM	4
217030	Cannabinoid CB <sub>1</sub>	401143	hum	2	10 µM	-15

Note: Items meeting criteria for significance (≥50% stimulation or inhibition) are highlighted.  
 \* Batch: Represents compounds tested concurrently in the same assay(s).  
 gp=Guinea pig, hum=Human

### Experimental Results


Cat #	Assay Name	Batch*	Spec.	Rep.	Conc.	% Inh.
217100	Cannabinoid CB <sub>2</sub>	401141	hum	2	10 µM	21
219500	Dopamine D <sub>1</sub>	401243	hum	2	10 µM	-14
219600	Dopamine D <sub>2L</sub>	401370	hum	2	10 µM	12
219700	Dopamine D <sub>2B</sub>	401241	hum	2	10 µM	13
219800	Dopamine D <sub>3</sub>	401243	hum	2	10 µM	9
224010	Endothelin ET <sub>A</sub>	401229	hum	2	10 µM	7
224110	Endothelin ET <sub>B</sub>	401366	hum	2	10 µM	-4
226010	Estrogen ERα	401171	hum	2	10 µM	75
226810	GABA <sub>A</sub> , Chloride Channel, TBOB	401172	rat	2	10 µM	-8
226600	GABA <sub>A</sub> , Flunitrazepam, Central	401246	rat	2	10 µM	-17
228510	GABA <sub>B</sub> , Non-Selective	401185	rat	2	10 µM	3
232030	Glucocorticoid	401208	hum	2	10 µM	4
232600	Glutamate, AMPA	401173	rat	2	10 µM	8
232700	Glutamate, Kainate	401296	rat	2	10 µM	-3
232810	Glutamate, NMDA, Agonism	401224	rat	2	10 µM	5
232910	Glutamate, NMDA, Glycine	401224	rat	2	10 µM	9
239300	Growth Hormone Secretagogue (GHS, Ghrelin)	401154	hum	2	10 µM	3
239610	Histamine H <sub>1</sub>	401249	hum	2	10 µM	6
239710	Histamine H <sub>2</sub>	401256	hum	2	10 µM	-8
239820	Histamine H <sub>3</sub>	401166	hum	2	10 µM	-4
243000	Insulin	401420	rat	2	10 µM	8
252200	Motilin	401156	hum	2	10 µM	13
252610	Muscarinic M <sub>1</sub>	401216	hum	2	10 µM	-1
252710	Muscarinic M <sub>2</sub>	401231	hum	2	10 µM	-2
252810	Muscarinic M <sub>3</sub>	401231	hum	2	10 µM	21
252910	Muscarinic M <sub>4</sub>	401217	hum	2	10 µM	14
258590	Nicotinic Acetylcholine	401227	hum	2	10 µM	-7
260130	Opiate δ <sub>1</sub> (OP1, DOP)	401150	hum	2	10 µM	-6
260210	Opiate κ (OP2, KOP)	401150	hum	2	10 µM	7
260410	Opiate μ (OP3, MOP)	401232	hum	2	10 µM	5
299005	Progesterone PR-B	401321	hum	2	10 µM	65
266700	Purinergic P2X	401151	rabbit	2	10 µM	3
268810	Purinergic P2Y	401152	rat	2	10 µM	-6
271110	Serotonin (5-Hydroxytryptamine) 5-HT <sub>1A</sub>	401182	hum	2	10 µM	10
271650	Serotonin (5-Hydroxytryptamine) 5-HT <sub>2A</sub>	401181	hum	2	10 µM	-2

Note: Items meeting criteria for significance (≥50% stimulation or inhibition) are highlighted.  
 \* Batch: Represents compounds tested concurrently in the same assay(s).  
 gp=Guinea pig, hum=Human


### Experimental Results

Cat #	Assay Name	Batch*	Spec.	Rep.	Conc.	% Inh.
271700	Serotonin (5-Hydroxytryptamine) 5-HT <sub>2B</sub>	401168	hum	2	10 µM	-4
271800	Serotonin (5-Hydroxytryptamine) 5-HT <sub>2C</sub>	401180	hum	2	10 µM	3
202000	Transporter, Adenosine	401137	gp	2	10 µM	-2
220320	Transporter, Dopamine (DAT)	401164	hum	2	10 µM	84
226400	Transporter, GABA	401186	rat	2	10 µM	-5
204410	Transporter, Norepinephrine (NET)	401164	hum	2	10 µM	16
274030	Transporter, Serotonin (5-Hydroxytryptamine) (SERT)	401182	hum	2	10 µM	3
287530	Vasopressin V <sub>1A</sub>	401219	hum	2	10 µM	-7


// IC<sub>50</sub>s were determined,  
 → Please refer to page 5




**Cpd 1a**  
Hydra  
WO 2015164643




**Cpd 2a**  
Glenmark  
WO 2013183035




**Cpd 3a**  
Eli Lilly  
WO 2019152465




**Cpd 4a**  
Roche  
WO 2018029288




**Cpd 5a**  
Abbvie  
\*Probe A-079




**Cpd 6a**  
Novartis  
J. Med. Chem., 2014, 57, 5129-5140




**Cpd 7a**  
AstraZeneca  
WO 2014184248




**Cpd 8a**  
Amgen  
J. Med. Chem., 2016, 59, 2794-2809




**Cpd 9a**  
Genentech Inc.  
WO 2018162607




**Cpd 10a**  
Almirall  
WO 2017060488




**Cpd 11a**  
Janssen  
WO 2014049047




**Cpd 12a**  
Orion  
WO 2014053694



**Cpd 13a**  
Merck Sharp & Dome  
WO 2011043954




**Cpd 14a**  
Ario Pharma  
WO2014135617




**Cpd 15a**  
Shanghai Jing'an Pharma.  
CN 108558831






Ca<sup>2+</sup>-transients of BAY-390



BAY-390 hTRPA1 IC<sub>50</sub> fit



- TRPA1 overexpressing cell line (e.g. CHO hTRPA1 GCaMP6)
- Add test compounds for 10 min
- Activation with agonist supercinnemaldehyde (SCA)
- Readout: calcium sensor GCaMP6 (or Fluo8)