



Science For A Better Life

# Donated Probe BAY-876


GLUT1 inhibitor

Charlotte Kopitz, Bernd Buchmann

# Introduction – Scientific Rationale

## Key Points

- Most cancers exhibit increased glucose uptake and aerobic glycolysis, that is, the Warburg effect
- Glucose transporters are rate-limiting checkpoints, abnormally regulated in cancer
- Metabolic rewiring and increased glycolysis drives resistance to SOC therapy
- Glucose transporter 1 (GLUT1 / SLC2A1) is a high-affinity glucose transporter with a defined expression pattern in normal tissue
- GLUT-1 is one of the major glycolysis players up-regulated by hypoxia-inducible factor 1 alpha (HIF1 $\alpha$ )




Carcinoma of the bladder, figure adapted from Reis et al. 2011



# Introduction – Scientific Rationale

## Rationale for a GLUT1 selective inhibitor


- Class I GLUT transporters have distinct physiological roles and share structural similarity
- GLUT-1 overexpression has been observed in multiple tumor indications
- GLUT-1 is one of the major glycolysis players up-regulated by hypoxia-inducible factor 1 alpha (HIF1 $\alpha$ )
- Other Glucose transporters are additionally expressed in heart, muscle, liver, kidney, pancreas (beta cells)
- panGLUT inhibitors have strong side effects in mice at efficacious dose (brain, heart)





Carcinoma of the bladder, figure adapted from Reis et al. 2011

**Specifically blocking GLUT1 should target tumors  
while minimizing the risk of side effects**

# Introduction – Rationale



Protein Class I	Km (mM)	Major sites of expression	Proposed function	Homology Human/ Mouse
GLUT1	3-7	Ubiquitous, overexpressed in tumors	Transport across blood tissue barriers	96.3
GLUT2	17	Liver, kidney, small intestine, $\beta$ -cells	„Glucose sensor“ in pancreatic $\beta$ -cells Transport across blood liver barrier	81.6
GLUT3	1.4	Neuronal cells	Neuronal transport	83.3
GLUT4	6.6	Muscle, fat, heart	Insulin-regulated transport	95.3




- class I and II: N-linked glycosylation site at 1st exofacial loop between TM1 and TM2
- class III: shorter loop 1; glycosylation site at larger loop 9

⇒ Selectivity against GLUT2, GLUT3 and GLUT4 regarded necessary for therapeutic window



# Probe BAY-876 – Profile



**BAY-876**

<i>Lead-like properties</i>	
MW <sub>corrected</sub>	441.2 g/mol
TPSA	127 Å <sup>2</sup>
Measured logD (pH 7.5)	2.7
Calculated logD (pH 7.5)	2.6
Solubility (pH 6.5)	0.8 mg/L
Solubility (Tween80/Soya 5:1)	1610 mg/L
Stability in r/h plasma	stable


<i>Highly potent, cellularly active and very selective</i>		
IC <sub>50</sub> GLUT-1	2 nM	
Selectivity	IC <sub>50</sub> GLUT-2	9.4 μM
	IC <sub>50</sub> GLUT-3	1.6 μM
	IC <sub>50</sub> GLUT-4	270 nM
IC <sub>50</sub> Cellular mechanistic assay (Glucose uptake inhib. / HeLa-MaTu)	3.2 nM	
Selectivity Kinases (# = 18)	> 20,000 nM	
Selectivity Eurofins Panel (# = 68)	> 5,000 nM Adenosine A3 : IC <sub>50</sub> = 1,140 nM Prostaglandin EP4 : IC <sub>50</sub> = 1,120 nM	

<i>Promising in vitro / in vivo pharmacokinetic profile</i>				
CL <sub>int</sub> Hepatocytes	Rat	1.9 L/h/kg	Rat PK in vivo	
	Mouse	0.93 L/h/kg	Cl <sub>b</sub>	0.33 L/h/kg
CL <sub>int</sub> Microsomes	Human	0.18 L/h/kg	V <sub>SS</sub>	0.79 L/kg
	P <sub>app</sub> AB	78 nm/sec	t <sub>1/2</sub>	4.9 h
Caco2	Efflux ratio	2.5 fold	F	85 %




# Potential Negative-Probe

## BAY-588 & BAY-480 – Basic Information



**BAY-588**



**BAY-480**

### Basic information

<b>IC<sub>50</sub> GLUT-1</b>		<b>1.2 μM</b>
<b>Selectivity</b>	<b>IC<sub>50</sub> GLUT-2</b>	<b>&gt; 10μM</b>
	<b>IC<sub>50</sub> GLUT-3</b>	<b>5.5 μM</b>
	<b>IC<sub>50</sub> GLUT-4</b>	<b>0.5 μM</b>
<b>MW<sub>corrected</sub></b>		<b>472.3 g/mol</b>
<b>TPSA</b>		<b>103 Å<sup>2</sup></b>
<b>Measured logD (pH 7.5)</b>		<b>3.9</b>
<b>Solubility (pH 6.5 / from DMSO sol.)</b>		<b>&lt; 1 mg/L</b>

### Basic information

<b>IC<sub>50</sub> GLUT-1</b>		<b>2.6 μM</b>
<b>Selectivity</b>	<b>IC<sub>50</sub> GLUT-2</b>	<b>&gt; 10μM</b>
	<b>IC<sub>50</sub> GLUT-3</b>	<b>&gt; 10 μM</b>
	<b>IC<sub>50</sub> GLUT-4</b>	<b>2.8 μM</b>
<b>MW<sub>corrected</sub></b>		<b>512.3 g/mol</b>
<b>TPSA</b>		<b>116 Å<sup>2</sup></b>
<b>Measured logD (pH 7.5)</b>		<b>2.4</b>
<b>Solubility (pH 6.5 / from DMSO sol.)</b>		<b>1.9 mg/L</b>




# Initial Pharmacology of BAY-876 in HeLa-MaTu Cells

## In vitro pharmacology of BAY-876


### Induction of cell death under hypoxia in 3D assay

HeLa MaTu, red= total spheroid, green= dead cells




Staurosporine control 100µM 10µM 1µM 100nM 10nM 1nM  
BAY-876

### Inhibition of glucose uptake




### Inhibition of proliferation (2D)



## BAY-876 in vivo

### Additive anti-tumor efficacy of BAY-876 in combination with radiation in vivo



- NMP/PEG 400 (1:9) 10 ml/kg QD p.o.
- ▲ vehicle 10 ml/kg QD p.o. + Radiation 10 Gy AL (4.2 min)
- BAY-876 2 mg/kg QD p.o.
- BAY-876 2 mg/kg QD p.o. + Radiation 10 Gy AL (4.2 min)

- \* P = 0.016 (t-test)
- \* P < 0.05 vs. vehicle (Dunnett's method)
- ↓ Radiation
- ↓ BAY-876 treatment



## Summary / Conclusion

- BAY-876 is a selective GLUT1 inhibitor meeting all chemical probe criteria:
  - Biochemical activity  $IC_{50}$  at 2 nM
  - On-target cellular activity (glucose uptake assay in HeLa-MaTu cell line) at 3.2 nM
  - Selectivity against family members was tested
- Selectivity in kinase and receptor panel was performed
- Additionally, two structure related compounds with low activity were identified and will be provided
- ***BAY-876 is a highly potent, cellular active and very selective inhibitor with a promising in vitro / in vivo pharmacokinetic profile***
- ***We would therefore provide the GLUT1 inhibitor BAY-876 and an inactive/negative control as a chemical probe***





# Acknowledgement

Andrea Hägebarth  
Andreas Bernthaler  
Arndt Schmitz  
Arwed Cleve  
Bernd Buchmann  
Carolyn Algire  
Charlotte Kopitz  
Eleni Lagkadinou  
Georg Kettschau  
Hanno Wild  
Hartmut Rehwinkel  
Heike Petrul  
Herbert Himmel  
Holger Siebeneicher

Ildikó Terebesi  
Iring Heisler  
Jens Geisler  
Judith Günther  
Karl Ziegelbauer  
Katja Prella  
Kirstin Meyer  
Luisella Toschi  
Marcus Bauser  
Marcus Koppitz  
Maria Quanz  
Mèlanie Hèroult  
Michael Brands  
Olaf Prien

Roland Neuhaus  
Roman Hillig  
Stephan Völkening  
Thomas Müller



Science For A Better Life

Thank you!