

## **Donated Chemical Probe**

## Chemical Probe BAY-707 MTH1 Inhibitor


March, 2018

Matyas Gorjanacz, Manuel Ellermann, Ashley Eheim \_\_\_\_\_





Scientific rationale: MTH1 as an anti-cancer target



// Cancer cells are characterized by oxidative
stress, which can damage dNTPs/DNA
// MTH1 (MutT homolog 1, NUDT1) prevents

oxidized dNTP incorporation into DNA by hydrolyzing 8-oxo-dGTP and 2-OH-dATP to 8oxo-dGMP and 2-OH-dAMP


// MTH1 expression and activity is up-regulated
in many cancers compared to normal tissue

// MTH1 is non-essential in normal cells (MTH1
knockout mice show only mild symptoms,
Tsuzuki, 2001)

// Initial disease hypothesis: inhibition of MTH1 will result in aberrant incorporation of oxidized nucleotides into DNA, subsequently leading to DNA damage, mutations, genomic instability and cancer cell death at excellent tolerability




*Initial supporting literature for MTH1* 



// MTH1 is overexpressed in cancer cells to oppose the DNA damage effects of high ROS levels // "Cancer-specific lethality" was described upon RNAi-mediated knockdown of MTH1 & upon treatment with small molecular weight MTH1 inhibitory tool compounds (e.g. TH588, (S)-crizotinib) // In vivo evidences also supported the assumption that MTH1 is required for cancer cell survival



## **Development of novel MTH1 inhibitors**




CL, clearance; IC50, half-maximal inhibitory concentration; logD, distribution coefficient; LM, liver microsomes; Sw, aqueous solubility.

Ellermann et al., 2017; ASC Chem Biol

- // Fragment-based screening and structure-based drug design led to discovery of a novel and selective MTH1 inhibitor with low nanomolar enzymatic activity ( $IC_{50}$ = 2.3+/-0.8 nM, n=6)
- // High selectivity in an in-house kinase panel, favourable physicochemical profile and promising *in vitro* pharmacokinetic properties with high metabolic stability and good cell permeability
- // Substrate competitive binding to the active site of MTH1







// Cellular Thermal Shift Assay (CETSA) used to demonstrate on-target cellular activity of BAY-707 and additional structurally related MTH1 inhibitors from the same compound class

// BAY-707 demonstrate a superior cellular target engagement (EC<sub>50</sub>= 7.6 nM) over the tool compound TH588 (EC<sub>50</sub>= 133 nM) // Good correlation between he biochemical potency and cellular target engagement of Bayer's MTH1 compound class // BAY-707 is a potent, selective and cellularly active MTH1 inhibitor with good PK properties; therefore it is suitable to validate the cellular functions of MTH1, e.g. the MTH1 cancer dependency




#### MTH1 is not required for cancer cell survival

Cell line <sup>a</sup>	Indication	TH588 Gl₅₀, μM <sup>b</sup>	BAY-707 GI <sub>50</sub> , μΜ <sup>b</sup>
HMEC	Normal breast	3.2	>30
NCI-H358	Lung cancer	4.9	>30
NCI-H460	Lung cancer	7.1	>30
A549	Lung cancer	4.0	>30
MCF7	Breast cancer	3.5	>30
MDA-MB-231	Breast cancer	6.3	>30
U2OS	Bone cancer	2.6	>30
HeLa	Cervical cancer	4.2	>30
SW480	Colon cancer	5.3	>30

aln total, more than 20 different cancer cell lines from varying indications and scientific rationales were tested.

 $^{b}$ Growth inhibition with the indicated compounds was performed in 6 day long assays. Both compounds were tested at concentrations of up to 30  $\mu$ M.

GI50, half-maximal growth inhibitory concentration; HMEC, human mammary epithelial cells.



// TH588 tool compound demonstrate equal cytotoxicity in normal and cancer cells

// BAY-707 demonstrate no cytotoxicity (neither in 2D nor 3D, independently of ROS levels, MTH1 expr.)


// No correlation between the biochemical potency and cytotoxicity of Bayer's MTH1 inhibitors

// No induction of double strand DNA breaks (DSBs) with BAY-707 and no correlation between γ-H2AX EC50 and biochemical potency of Bayer's MTH1 inhibitors

// DSBs observed with some MTH1 inhibitors are independent of their enzymatic activity and likely due to off-target effects



MTH1 is not essential for sanitization of oxidized dNTPs



#### DNA incorporation of 8-oxo-dGTP in U2OS cells

Ellermann et al., 2017; ASC Chem Biol


// Genomic incorporation of 8-oxo-dGTP was measure in HCA of immunostained samples
// BAY-707 and structurally related MTH1 inhibitors do not result in increased nuclear incorporation of
 damaged nucleotides (measured in γ-H2AX and 8-oxo-dGTP assays)
// Off target substativisity of TUE 88 and MTU1 siDNAs are the primary reason for the increased pusclear


// Off-target cytotoxicity of TH588 and MTH1 siRNAs are the primary reason for the increased nuclear incorporation of 8-oxo-dGTP

// In living cells MTH1 is not essential for sanitization of oxidized nucleotides



MTH1 inhibition in mono- or combination-therapies has no effect on in vivo tumor growth





Ellermann et al., 2017; ASC Chem Biol

- // No significant body weight loss observed with BAY-707 and other Bayer MTH1 inhibitors
- // Exposure for both BAY-707 covered biochemical  $IC_{50}$  and CETSA  $EC_{50}$  for ~24h
- // MTH1 inhibition in mono- or combination-therapies has no effect on in vivo tumor growth



Compound Comparison to MTH1 Negative Control BAY-604

Property	BAY-707	BAY-604	TH588
IC <sub>50</sub> MTH1 [nM]	2.3	>20000	12
logD@ pH 7.5	1.6	1.9	2.5
Sw pH 6.5 [mg/L]	288	2170	38
Caco-2 A-B [nm/s] (efflux ratio)	148 (1.4)	68 (3.2)	307 (0.5)
Rat hep. CL [L/h/kg] (Fmax)	0.54 (87%)	0.41 (90%)	1.92 (54%)
Hum. LM CL [L/h/kg] (Fmax)	0.29 (78%)	0.015 (99%)	0.64 (51%)
GI <sub>50</sub> [μΜ] SW480	>30	>30	5.3
GI <sub>50</sub> [µM] NCI-H358	>30	>30	4.8

// BAY-707 is a low nanomolar MTH1 inhibitor with  $IC_{50}$  = 2.3+/-0.8 nM enzymatic activity

// BAY-604 is structurally related compound to BAY-707 with similar PhysChem and PK properties

// BAY-604 is the negative control of BAY-707 and demonstrates no enzymatic activity until the highest does tested (20 μM)

// BAY-604 demonstrate no cancer cell growth inhibition up to the highest concentration tested (30  $\mu$ M)



// BAY-707 is an MTH1 inhibitor fulfilling all criteria for a chemical probe:

// Low nanomolar biochemical potency (IC<sub>50</sub>= 2.3 nM)

- // Substrate-competitive binding
- // Good membrane permeability and single-digit cellular target engagement (CETSA assay EC<sub>50</sub>= 7.6 nM)
- // Selective against an in-house kinase panel. In contrast to the currently widely used MTH1 tool compound TH588, BAY-707 demonstrate no off-target-related cytotoxicity

// BAY-707 de-validated MTH1 as a broad-spectrum non-oncogenic cancer dependency

BAY-707 will allow to further study the biology of MTH1 in cell cultures and living organisms without a limitation of off-target-related cytotoxicity

BAY-604 is a negative control of BAY-707

We ask for acceptance of MTH1 inhibitor BAY-707 as donated chemical probe





Matyas Gorjanacz\* Manuel Ellermann\* Ashley Eheim\* Maria Quanz Andrea Glasauer Anja Richter Stefanie Bunse Katrin Nowak-Reppel Horst Irlbacher Benjamin Bader Judith Guenther Hanna Meyer Roland Neuhaus Jörg Weiske

**Sprint Bioscience** Fredrik Ranm<sup>\*</sup> Jenny Viklund Tobias Ginman Fredrik Rahm\* Johan Lindström Martin Andersson Andrea Haegebarth Ingo Hartung Marcus Bauser Anke Mueller-Fahrnow Ursula Egner

Jan Huebner

\* Core Team



# Thank You





### Further cellular validation of MTH1

