



Donated Chemical Probe

Chemical Probe BAY-826
Tie/DDR Inhibitor

March, 2018

Sylvia Gruenewald, Lars Baerfacker





Tie/DDR probe BAY-826:

Scientific rationale: Tie2 as an anti-cancer target

RTK Tie2 as a “classical” anti-angiogenesis target

- // Ang1/2-Tie2 signaling: important role in angiogenesis and vessel maturation
- // Tie2 inhibition: impairs angiogenesis & reduces tumor growth
 - shows combination benefit with anti-VEGFR therapy

Tie2 as anti-tumor cell target beyond angiogenesis

- // Survival of **Tie2-positive AML cells** sustained through autocrine Ang1/Tie2 loop
- // **Tie2-positive hematological tumor cells** may adhere to the bone marrow niche thus being protected from chemotherapy
- // Activation of **gliomal Tie2** increases tumorigenesis and invasive phenotype *in vivo*
- // **Tie2** activation in **gliomal and brain tumor stem cells** contributes to chemoresistance
- // Infiltration of **Tie2-expressing macrophages** (TEMs) is implicated to promote angiogenesis and **metastatic dissemination**, potentially of relevance in various tumor indications incl. HCC.




Tie/DDR probe BAY-826:

Disease hypothesis based on literature data


Angiogenesis - Tie2 function is essential for development of embryonic vasculature, for maintenance of blood vessels and for the maturation of newly formed vessels

In situ BS-1 lectin fluorescent staining of corneal limbus vessels 6 days after pellet insertion



Asahara, T. *et al.* (1998) *Circ Res* 83, 233-240

Comparison of the sTie2 secreting A375 melanoma xenografts with parental cells with and without treatment of PTK787/ZK222584



Siemeister *et al.* (1999) *Canc. Res.* 59, 3185 and WO2001/097850


// sTie2-Fc precluded modulation of VEGF-induced neovascularization by either angiopoietin (Ang2 as well as Ang1) (left)
 // Tie2 inhibition on top of VEGFR2 inhibition improves efficacy




Tie/DDR probe BAY-826:

Disease hypothesis based on literature data

AML - Tie2 & Ang1 expressed in subset of AML blasts



Blocking of Ang1-Tie2 signaling in Tie2+ but not Tie2- AML cells increases apoptosis




a) Riccioni *et al.* (2007) *Stem Cells*, 25, 1862. b) Arai *et al.* (2004) *Cell*, 118, 149. c) Ichihara *et al.* (2011) *BBRC*, 416, 239.

Glioma - Tie2 expression in neoplastic glial cells correlates with grade


Tie2 staining on human TMA sections

normal brain (NB), low grade astrocytoma (LGA), anaplastic astrocytoma (AA), and glioblastoma multiforme (GBM)



Activation of gliomal Tie2 increases tumorigenesis and invasive phenotype *in vivo*

U251.vector or U251.Tie2 cells were injected with ECs into the brains of immunocompromised mice → infiltrative/multifocal component of tumors in Tie2+/EC




a) Lee *et al.* (2006) *Mol. Canc. Res.*, 4, 915. b) Liu *et al.* (2010) *Oncotarget*, 1, 700. c) Martin *et al.* (2009) *Oncogene*, 28, 2358.



Tie/DDR probe BAY-826:

Overall profile



▪ Molecular Properties

MW [g/mol]	559
MWcorr [g/mol]	490
TPSA [Å ²]	88
Rotatable bonds	5

▪ PhysChem

Sw ^{pH 6.5} [mg/L]	1.8
log D (pH 7.5)	3.61

▪ Pharmacology

Tie2 K _D (KINOMEScan™)	1.6 nM
Tie2 IC ₅₀ (in-house kinase assay)	0.45 nM
Tie2 IC ₅₀ (HUVEC pTie2-ELISA)	1.3 nM
Tie1 K _D (KINOMEScan™)	0.9 nM
DDR1 K _D (KINOMEScan™)	0.4 nM
DDR2 K _D (KINOMEScan™)	1.3 nM
VEGFR2 K _D (KINOMEScan™)	1.6 μM
VEGFR3 IC ₅₀ (in-house kinase assay)	0.44 μM
FGFR1/3 IC ₅₀ (in-house kinase assay)	>10 μM
PDGFR-β IC ₅₀ (in-house kinase assay)	>20 μM

▪ In vitro PK

		Clint [L/h/kg]	Fmax [%]
LM	Human	0.43	68
	Mice		
	Rat	1.91	54
Hep	Rat	0.96	77
CaCo2	A-B [nm/s]	78	Ratio
	B-A [nm/s]	55	0.7


▪ Safety (LeadProfilingScreen, total # of assays:68)

BAY-826 (10 μM)	Adenosine A3 (65 % inh.), Opiate κ (75% inh.), hERG (66 % inh.), Sodium Channel site2 (84% inh.)
---------------------------	--



Tie/DDR probe BAY-826:

Negative control *BAY-309*



Pharmacology

Tie2 K_D	>20 μM
------------	-------------

Molecular Properties

MW [g/mol]	519
MWcorr [g/mol]	450
TPSA [\AA^2]	64
Rotatable bonds	5

PhysChem

$S_w^{pH\ 6.5}$ [mg/L]	0.14
log D (pH 7.5)	3.71

In vitro PK

CaCo2	A-B [nm/s]	B-A [nm/s]	Ratio
	No data (below detection limit)		



Tie/DDR probe BAY-826:

Biochemical and cellular potency and selectivity

Biochemical activity*	K _D [nM]
Tie2	1.6 (0.5, IC ₅₀)
Tie1	0.9
DDR1	0.4
DDR2	1.3
LOK	5.9
EPHB6	25
LYN	40
MERTK	66 (IC ₅₀)
ABL1-non-phosphorylated	38
ABL1-phosphorylated	510
ABL1 (T315I)-non-phosphorylated	100
ABL1 (T315I)-phosphorylated	310
KIT	150 (>20000, IC ₅₀)
BRAF (V600E)	370
BRAF	730
RAF1	1400
VEGFR2	1600 (118, IC ₅₀)
VEGFR3	439 (IC ₅₀)
FGFR1	16200 (IC ₅₀)
FGFR3	12800 (IC ₅₀)
PDGFR-β	>20000 (IC ₅₀)

*K_D values determined by DiscoverX Corp. on the KINOMEScan platform

Summary of cellular NanoBRET™ assay data @SGC

	BAY-826 IC ₅₀ (normed to Tie-2 data)*	Ratio of IC ₅₀ s (BAY-309/BAY-826)*
Tie2	1.0	1604.7
DDR2	4.5	3056.1
DDR1	1.4	369.9
Tie1	8.1	23.4
STK10***	317.0	
EPHB6	5822.8	2.9

* Dose response curves and derived IC₅₀ values, see backup section

// BAY-826 binds to Tie1, Tie2, DDR1 and DDR2 with a K_D of ~ 1 nM

// BAY-826 is selective versus the known angiogenic RTKs VEGFR1/2/3, FGFR1/2/3/4 and PDGFR-α/β

// BAY-826 is equipotent vs DDR1/2 (and Tie1 kinase) in cellular assays and more than three orders of magnitude less active vs STK10 and EPHB6 kinase

// The neg. control BAY-309 reveals a sufficient cellular selectivity vs BAY-826 for all six kinases tested



Tie/DDR probe BAY-826:

Summary / conclusion

Probe criteria	
Inhibitor/agonist potency: goal is < 100 nM (IC ₅₀ , Kd)	Surpasses criteria; biochemical assay (Tie-2) with IC ₅₀ 0.45 nM;
Selectivity within target family: goal is > 30-fold	Surpasses criteria; selectivity > 100 fold vs all other angiogenic kinases, (# of kinases tested 453) albeit cellular equipotency vs DDR1/2 and Tie-1
Selectivity outside target family: describe the off-targets (which may include both binding and functional data)	Surpasses criteria; Promising LeadProfilingScreen: 4 out of 68 assays hit (all in μM potency range)
On target cell activity for cell-based targets: goal is < 1 micromolar IC ₅₀ /EC ₅₀	Surpasses criteria; functional cellular assays: pTie2-ELISA, HUVEC-cells with IC ₅₀ 1.3 nM; NanoBRET™ assay with IC ₅₀ 0.7 nM
On target cell activity for secreted targets: appropriate alternative such as mouse model or other mechanistic biological assay, e.g., explant culture	n/a
Neg ctrl: <i>in vitro</i> potency - > 100 times less; Cell activity - > 100 times less potent than the probe	Surpasses criteria; > 10,000 times less in biochemical assay (Tie-2); > 1,000 times less in cellular assay (DDR1/2, STK10, EPHB6)

We ask for acceptance of Tie/DDR inhibitor BAY-826 as chemical probe, accompanied by BAY-309 as negative control which may allow to selectively study the biology of Tie/DDR signaling *in vitro* and *in vivo* as it does not target other angiogenic kinases



Tie/DDR probe BAY-826:

Acknowledgement

Sylvia Grünewald
Lars Bärfacker
Stefan Kaulfuss
Dieter Heldmann
Mario Lobell
Michael Härter
Frank Süßmeier
Ulf Bömer
Martina Schäfer
Jörg Weiske
Wolfgang Halfbrodt
Bernd Buchmann
Ingo Hartung
Alexander Walter
Anna-Lena Frisk

Karl Ziegelbauer
Andrea Haegebarth
Marcus Bauser
Rolf Jautelat
Holger Hess-Stumpp
Anke Mueller-Fahrnow
Ursula Egener

Cora Scholten

Oncotest

Julia Schüler
Kerstin Klingner



Thank You





Tie/DDR probe BAY-826 & negative control BAY-309:

NanoBRET™ assay data, SGC results

170829/170901 NanoBRET Tie2 Probe and Neg Ctrl

