




Donated Chemical Probe

*Complex 1 Inhibitor
Probe BAY-179*

November, 2019

Sven Christian, Jeffrey Mowat, Roman Hillig






Complex 1 probe BAY-179

Target rationale and cancer


Physiological role of complex-1



- The electron transport chain is located in the inner mitochondrial membrane
- It transfers electrons from electron donors to electron acceptors via redox reactions (Complex I-IV), which is coupled to the transfer of protons (H⁺ ions) across the membrane. This creates an electrochemical proton gradient that drives ATP synthesis (Complex V). The final acceptor of electrons in the electron transport chain is molecular oxygen.
- Mitochondria are key regulators of energy supply & cell death

Target Rationale

- Oncogenic mutations have been described to sensitize tumor cells to Complex I inhibition (e.g. LKB1 -/-)



- Inhibition of BRAF^{V600E} by Vemurafenib increases PGC1a, MITF and OXPHOS [Haq et al., 2013, Cancer Cell] and therefore sensitizes to Complex I inhibition
- A subset of diffuse large B-cell lymphoma (Oxphos-DLBCL) has been characterized by enhanced Oxphos [Caro et al., 2012, Cancer Cell]
- Former HIF-FR project showed that Complex 1 inhibitors have enhanced efficacy when used in combination with Radiation

Genetic mutations and consequent metabolic changes have been shown to sensitive cancer cells to Complex 1 inhibition




Complex 1 probe BAY-179

Known Complex-1 inhibitors


Metformin

- Oral antidiabetic drug in the biguanide class. The first-line drug of choice for the treatment of type 2 diabetes
- approved in Canada in 1972, approval by FDA for type 2 diabetes in 1994
- $IC_{50} > 10 \mu M$
- slow membrane-potential-driven accumulation of the drug within the mitochondria




Phenformin

- Biguanide class
- marketed as DBI by Ciba-Geigy, but was withdrawn from most markets in the late 1970s due to a high risk of lactic acidosis, which was fatal in 50% of cases
- $IC_{50} > 1 \mu M$



Rotenone

- broad-spectrum insecticide, piscicide, and pesticide
- $IC_{50} = 20 \text{ nM}$




- Rotenone is neurotoxic and induces a Parkinson-like phenotype in mice
- In house data show that toxicity can not be rescued by expression of NDI1 => Complex 1 independent effect

No potent and selective Complex 1 inhibitor commercially available



Complex 1 probe BAY-179

Overall profile



Pharmacology

IC₅₀ H1299 HTS ATP (human)	79 nM
IC₅₀ H1299 CTG ATP (human)	33 nM
IC₅₀ 4T1 CTG ATP (mouse)	38 nM

Safety

hERG IC₅₀ volt clamp	3.0 μM
Off-targets	Serotonin 5HT2b in-house data Agonism: no effect
	Antagonism: IC₅₀ = 600 nM (sensitive set-up, 10nM 5HT agonist)

Molecular Properties

MWcorr [g/mol]	416
i.s.o. PS (0-10)	3
TPSA [Å²]	71
clogP / cLogD@pH7.5	4.0/3.6
Rotatable bonds	4

PhysChem

Sw solid @pH 6.5 [mg/L]	1.0
Log D (pH 7.5)	2.8
Stability (plasma h,r,m)	stable
Stability (pH 1,7,10) [% remaining after (1h)]	100/100/100


In vitro PK

		CL_B [L/h/kg]		Fmax [%]			
LM	Human	0.55		58			
	Mouse (f)	1.9		65			
	Rat	0.77		82			
	Dog	0.53		75			
Hep	Rat	0.81		81			
	Dog	1.1		47			
	Mouse (f)	1.5		72			
	Human	0.64		52			
CaCo2		A-B [nm/s]		B-A [nm/s]		Ratio	
		92		41		0.44	
Cyp Inh profile [μM]		1A2	2C8	2C9	2D6	3A4	3A4 preinc
		> 20	7.8	15	> 20	> 20	> 20
Cyp Ind [μM]		PXR EC₅₀		PXR Rating		CYP Ind NOEL	
		>2μM (7.8% Emax)		green		3A4: NOEL 123μg/l loss of cell viability @ higher doses	



Complex 1 probe BAY-179

Cellular data - proliferation



Indication	Cell line	Patient-derived	IC ₅₀ [µM]	IC ₇₀ [µM]
Bladder	BXF135 2	x	0,023	0,048
Gastric	MKN1		0,034	0,085
Bladder	BXF121 8	x	0,065	0,104
Colon	LS174T		0,043	0,109
Gastric	GXF25 1	x	0,081	0,147
Lung	LXFA62 9	x	0,066	0,195
Ovary	BG1		0,142	0,207
Colon	Caco-2		0,129	0,209
Melanoma	SK- MEL-28		0,075	0,215
Colon	RKO		0,153	0,251
Melanoma	LOX		0,126	0,335

Unbiased approach to identify more sensitive tumor models (oncotest):

- Several cancer cell lines were identified to be highly sensitive towards Complex 1 inhibition
- Selected cell lines were tested in house in vivo
- Bioinformatic analysis identifies KRASmut to be enriched in resistant tumor cell lines



Complex 1 probe BAY-179

DMPK data of BAY-179


▪ *In vitro* PK

		LM CL _{blood} [L/h/kg]		HEP CL _{blood} [L/h/kg]		
Metabolic Stability	Mouse	1.9		1.5		
	Rat	0.77		0.81		
	Dog	0.53		1.1		
	Human	0.55		0.64		
Fraction unbound [%]	Williams_E	Mouse	Rat	Dog	Human	
	24	0.43	0.26	0.37	0.48	
Caco-2	A-B [nm/s]	B-A [nm/s]		Ratio		
	92	41		0.44		
CYP Inh [μM]	1A2	2C8	2C9	2D6	3A4	3A4 preinc
	>20	7.8	15	>20	>20	>20
PXR Assay	Rating:		green			
CYP Induction	1A2: NOEL 41μg/l		3A4: NOEL 123μg/l			loss of cell viability @ higher doses

- *In vitro* CL: low to moderate in hepatocytes and liver microsomes
- High protein binding in all species
- *In vivo* rat PK: Low CL_{blood}, high V_{ss}, long t_{1/2}, high oral bioavailability

▪ *In vivo* PK

Species		Rat	
Route		iv inf.	po
Dose	[mg/kg]	0.2	0.6
AUC _{norm}	[kg·h/L]	4.1	3.1
Cmax _{norm}	[kg/L]		0.23
t _{max}	[h]		4 to 7
CL _{matrix}	[L/h/kg]	0.25	
CL _{blood}	[L/h/kg]	0.40	
V _{ss}	[L/kg]	2.2	
t _{1/2}	[h]	6.2	6.2
F	[%]		76



BAY-179 is suitable for in vivo studies



Complex 1 probe BAY-179

GPCR data of BAY-179

The results of the current Eurofins-Cerep GPCR-Panel:

BAY-179

A2B (h) (antagonist effect) 57.8%

H2(h) (antagonist effect) 64.4%

MC4(h) (antagonist effect) 67.7%



M1(h) (antagonist effect) 53.6%

μ (MOP) (h) (antagonist effect) 51.3%

5-HT2B(h) (antagonist effect) 73.5%

BAY-070

5-HT2B(h) (antagonist effect) 59.9%






Complex 1 probe BAY-179


Bacterial Complex 1 as crystallization surrogate (Prof. Leonid Sazanov, IST Austria)

Structure of the entire *Thermus* complex I




Baradaran, Berrisford, Minhas and Sazanov, *Nature*, 2013

Probe candidate: **BAY-179**




Chemically related **compound-1**
(structure cannot be disclosed)

Dose response in biochemical activity assay with bacterial complex I*:



IC50 (human) = 79 nM
IC50 (bact.) = 5.3 μM
→ too weak for soaking



IC50 (human) = 5.1 nM
IC50 (bact.) = 192 nM
→ Sufficiently potent for soaking

- In collaboration with Prof. Leonid Sazanov, BAY-179 and related compounds were tested for activity on Complex 1 from *Th. Thermophilus* in a biochemical activity assay*.
- Some compounds showed very similar IC50 values in human vs. bact. assay, others showed differences (more often weaker on bact. than human enzyme).
- Probe BAY-179 itself showed poor potency, but chemically related compound-1 showed acceptable potency in bacterial surrogate system and was successfully soaked into crystals of bacterial Complex 1.


Crystal structure shows density for ligand in quinone-binding site (structure cannot yet be disclosed, manuscript in preparation).

* Bacterial activity assay: NADH:decylubiquinone assay; rate of NADH oxidation by complex I was followed at different compound concentrations which was subsequently used to plot a dose response curve



Complex 1 negative control BAY-070

Overall profile



Pharmacology

IC ₅₀ H1299 HTS ATP (human)	>30 μM
IC ₅₀ H1299 CTG ATP (human)	
IC ₅₀ 4T1 CTG ATP (mouse)	

Kinase Panel

Bayer Kinase Panel 34 Kinases Tested	All > 20 μM
---	-------------

Molecular Properties

MWcorr [g/mol]	429.5
i.s.o. PS (0-10)	2
TPSA [Å ²]	76
clogP / cLogD@pH7.5	2.6 / 2.3
Rotatable bonds	4

PhysChem

Sw solid @pH 6.5 [mg/L]	
Log D (pH 7.5)	2.1
Stability (plasma h,r,m)	
Stability (pH 1,7,10) [% remaining after (1h)]	

In vitro PK

		CL _B [L/h/kg]			Fmax [%]	
LM	Human					
	Mouse (f)					
	Rat					
Hep	Rat					
	Dog					
	Human					
CaCo2	A-B [nm/s]	127		B-A [nm/s]	245	
					Ratio 1.9	
Cyp Inh profile [μM]	1A2	2C8	2C9	2D6	3A4	3A4 preinc
Cyp Ind [μM]	PXR EC ₅₀	PXR Rating		CYP Ind NOEL		



Complex 1 probe BAY-179

Summary / conclusion

Probe criteria	
Inhibitor/agonist potency: goal is < 100 nM (IC ₅₀ , Kd)	not applicable – see on target cell activity
Selectivity within target family: goal is >30-fold	Surpasses criteria ; compound class shows selectivity against other complexes of the respiratory chain
Selectivity outside target family: describe the off-targets (which may include both binding and functional data)	Surpasses criteria ; clean LeadProfilingScreen; selectivity > 500 fold against numerous kinases tested
On target cell activity for cell-based targets: goal is < 1 micromolar IC₅₀/EC₅₀	Surpasses criteria ; Active in cellular mechanistic assay demonstrating cellular target engagement cellular proliferation assays; (IC ₅₀ = 79nM)
On target cell activity for secreted targets: appropriate alternative such as mouse model or other mechanistic biological assay, e.g., explant culture	Surpasses criteria ; suitable pharmacokinetic profile for in vivo studies
Neg ctrl: in vitro potency – > 100 times less; Cell activity – >100 times less potent than the probe	Surpasses criteria ; BAY-070 inactive (IC ₅₀ >20μM) in cellular assay

We ask for acceptance of Complex 1 inhibitor BAY-179 as chemical probe, accompanied by BAY-070 as negative control



Complex 1 probe BAY-179

Acknowledgement

MedChem Berlin	Hue-Quan Lao Andrea Koehlke-Nork Michael Tietz Anni Knitter Lisa Schällicke Martin Schötz Karl Sauvageot-Witzku	DMPK	Michaela Bairlein Karsten Denner Roland Neuhaus Michael Niehues
		BIP	Wolfgang Halfbrodt Stephan Völkening Jochen Kossmann
MedChem Wuppertal	Carl Nising Nils Griebenow	FTO Analysis	Peter Spreyer
MedChem external	Peakdale Support	Toxicology	Anna-Lena Frisk Ute Bach
TRG Oncology	Laura Schöckel Melanie Appel Maria Spelling Seren Nesan Michael Reinhardt Elisabeth Krahl Ulrike Lenhard	Exp. Med.	Eleni Lagkadinou
		ER moderators	Andrea Haegebarth Marcus Bauser [Holger Hess-Stumpp]
LD	Nina Scheweling Patrick Steigemann Roman Hillig	Structural Biology	Leonid Sazanov, IST Vienna, Austria






Thank You





Complex 1 probe BAY-179

In vivo studies in LKB1^{-/-} tumor model NCI-H460




Unfortunately, BAY-179 shows no *in vivo* efficacy was reported in LKB1^{-/-} tumor model NCI-H460



Complex 1 probe BAY-179


Synthesis scheme





Complex 1 negative control BAY-070

Synthesis scheme





Complex 1 probe BAY-179

Selectivity against other Complexes

Compound		Related Analog to BAY-179	Rotenone
Cellular ATP assay*	Complex I	25 nM	10 nM
Bio-chemical bovine OxPhos activity assay**	Complex I	47 nM	9,2 nM
	Complex II	> 30 µM	> 30 µM
	Complex III	> 30 µM	> 30 µM
	Complex IV	> 30 µM	> 30 µM
	Complex V	> 30 µM	> 30 µM

Analog of BAY-179 displays target specificity in biochemical assays

* Measured by cell based primary HTS assay (human cell line H1299 with ter luc ATP sensor +/- succinate)

** Measured by Abcam MitoTox complete oxphos activity assay kit using bovine isolated mitochondria



Complex 1 probe BAY-179

Complex 1 cellular activity assay description

Cellular Complex I activity

For the measurement of cellular Complex I activity, cellular ATP levels in absence or presence of the Complex II substrate succinate are measured. Therefore, H1299 cells (ATCC CRL-5803) were stably co-transfected with a pcDNA3 vector (purchased from Invitrogen™) encoding for ATP dependent Pyrearinus termitillumins larval click 20 beetle luciferase under control of a CMV promoter (H1299 tluc cell). In brief, H1299tluc cells (4000/well) are seeded into white 384 well plates. After one day in culturing in DMEM (2% FCS) without glucose, but supplemented with 11 mM galactose, 10 µL of a luciferin/test compound mixture (150 µM D-luciferin, 0.4% DMSO final concentration in 2 mM Ca²⁺ Tyrode +0.01% BSA) is added to each well and incubated for 1 h at 37° C, 25 5% CO₂ and luminescence measurement is performed. After this measurement 20 µL succinate (0.67 M, pH 5.3 in Tyrode, final concentration 25 mM) is added. The plate is then incubated for another 1 h at room temperature before the second measurement is performed.

The succinate dependency of cellular ATP level obtained in the second measurement gives information on Complex I specificity of compound action (i. e. Complex I inhibitors will lose their effect).

Measurements were always performed in quadruple. For determination of inhibition constants IC₅₀ and the Hill slope, the normalized luminescence intensities (DMSO control = 100%) were plotted as a function of the concentration of test compound, and data points were fitted by nonlinear regression with the equation

$$Y = \text{Bottom} + (\text{Top} - \text{Bottom}) / (1 + 10^{((\log \text{IC}_{50} - X) * \text{HillSlope})})$$

where Y is the measured luminescence normalized to DMSO control (=100%), X is the log of compound concentration and Top/Bottom are the plateaus of the curve (normalized luminescence units).